



# ***PRODUCT RANGE***

## **SPECIFICATION SHEETS**

March 2019

**CONTACT**



Telephone: +44 (0) 1484 865586, Email: [info@titanmed.co.uk](mailto:info@titanmed.co.uk), Web: [www.titanmed.co.uk](http://www.titanmed.co.uk)  
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW  
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288

# PERFORMANCE

## 2 & 3 SECTION



Our robust entry level Performance range of treatment couches are available in either 2 or 3 sections with electric height operation, making them ideal for a wide range of healthcare environments and treatments procedures.

**2 Section: TMP2E** - Electric height elevation

**3 Section: TMP3E** - Electric height elevation

### STANDARD SPECIFICATION

- Safe working load and lifting capacity of 170kg (26 stone)
- Excellent height range from 50cm (20") to 95cm (37")
- Hoist clearance along entire length of couch
- 125mm diameter individually locking castors
- 50mm high density flame-retardant foam
- Antimicrobial upholstery
- Fully compliant and CE marked
- Hand switch fitted as standard
- Gas spring assisted back rest from -30° to +90°

### OPTIONAL ACCESSORIES

- Single foot switch operation
- Paper towel holder
- Phlebotomy arm
- Breathing hole

### DIMENSIONS

- Length: 190cm
- Width: 66cm
- Safe Load: 170kg
- Height range: 50cm-95cm

### WARRANTY

- 1 year warranty on moving components
- 3 year warranty on frame

### STANDARDS


- EEC classification - MDD93/42/EEC Annex IX Class I
- Electrical standard - IEC601 Part 1 + 60601. P (EN)
- IP rating - IPx4
- Vinyl flame retardant to CRIB5 (BS7176 inc. foam)
- Foam to BS 5852: Part 2: 1982, Ignition Source 5 (CRIB5)



Safe working load and lifting capacity of 170kg (26 stone)




#### Performance 2 Section:



126cm

64cm



Height Range 50cm - 95cm


#### Performance 3 Section:



68cm

50cm

69cm



Height Range 50cm - 95cm



Telephone: +44 (0) 1484 865586, Email: [info@titanmed.co.uk](mailto:info@titanmed.co.uk), Web: [www.titanmed.co.uk](http://www.titanmed.co.uk)  
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW  
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288

# 2 SECTION ELECTRIC & MANUAL



A very strong, robust and aesthetically pleasing treatment couch available with optional breathing hole and removable plug. Our two section couches feature 55mm retractable positioning castors with an adjustable self-levelling foot making them ideal for a wide range of healthcare environments and treatments procedures.

**TMU2H** - Hydraulic height elevation    **TMU2E** - Electric height elevation

## STANDARD SPECIFICATION

- Safe working load and lifting capacity of 250kg (39 stone)
- Choice of 2 models - hydraulic or electric height elevation
- Excellent height range from 48cm to 97cm
- 50mm high density flame-retardant foam
- Antimicrobial contoured upholstery
- Fully compliant and CE marked
- The power assisted head section incorporates anti-trap design
- Hand switch fitted as standard on electric models
- Hoist clearance along entire length of couch
- Purpose made precision bearings on all pivot points
- Gas spring assisted back rest from -30° to +85°

Safe working  
load and lifting  
capacity of 250kg  
(39 stone)



## OPTIONAL ACCESSORIES

- Motorised backrest
- Single foot switch operation for electric models
- Double foot switch operation for electric models
- Paper towel holder
- Cot sides (powder coated, chrome plated or padded)

## DIMENSIONS

- Length: 192cm (head 69cm, seat 122cm)
- Width: 65cm
- Safe Load: 250kg
- Height range: 48cm - 97cm

## WARRANTY

- 3 year warranty on moving components
- Lifetime warranty on frame


## STANDARDS

- EEC classification - MDD93/42/EEC Annex IX Class I
- Electrical standard - IEC601 Part 1 + 60601. P (EN)
- IP rating - IPx4
- Vinyl flame retardant to CRIB5 (BS7176 inc. foam)
- Foam to BS 5852: Part 2: 1982, Ignition Source 5 (CRIB5)



122cm

69cm



Height Range 48cm - 97cm



Telephone: +44 (0) 1484 865586, Email: [info@titanmed.co.uk](mailto:info@titanmed.co.uk), Web: [www.titanmed.co.uk](http://www.titanmed.co.uk)  
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW  
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288

# 3 SECTION ELECTRIC & MANUAL



A very strong, robust and aesthetically pleasing treatment couch available with optional breathing hole and removable plug. Our three section plinths feature 55mm retractable positioning castors with an adjustable self-levelling foot making them ideal for a wide range of healthcare environments and treatments procedures.

**TMU3H** - Hydraulic height elevation

**TMU3E** - Electric height elevation

## STANDARD SPECIFICATION

- Safe working load and lifting capacity of 250kg (39 stone)
- Choice of 2 models - hydraulic or electric height elevation
- Gas spring assisted back rest from -30° to +85°
- The power assisted head section incorporates anti-trap design
- Excellent height range from 48cm to 97cm
- Purpose made precision bearings on all pivot points
- Hoist clearance along entire length of couch
- Hand switch fitted as standard on electric models
- 50mm high density flame-retardant foam
- Antimicrobial contoured upholstery
- Fully compliant and CE marked

## OPTIONAL ACCESSORIES

- Motorised backrest
- Single foot switch operation for electric models
- Double foot switch operation for electric models
- Paper towel holder
- Cot sides (powder coated, chrome plated or padded)

## DIMENSIONS

- Length: 192cm (head 69cm, seat 50cm, foot 68cm)
- Width: 65cm
- Safe Load: 250kg
- Height range: 48cm - 97cm


## WARRANTY

- 3 year warranty on moving components
- Lifetime warranty on frame

## STANDARDS

- EEC classification - MDD93/42/EEC Annex IX Class I
- IP rating - IPx4
- Electrical standard - IEC601 Part 1 + 60601. P (EN)
- Vinyl flame retardant to CRIB5 (BS7176 inc. foam)
- Foam to BS 5852: Part 2: 1982, Ignition Source 5 (CRIB5)

Safe working  
load and lifting  
capacity of 250kg  
(39 stone)



Telephone: +44 (0) 1484 865586, Email: [info@titanmed.co.uk](mailto:info@titanmed.co.uk), Web: [www.titanmed.co.uk](http://www.titanmed.co.uk)  
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW  
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288

# NEUROLOGY COUCH



A very strong, robust and aesthetically pleasing treatment couch with breathing hole and removable plug fitted as standard. Our two section neurology couch features 55mm retractable positioning castors with locking function and an adjustable self-levelling foot making them ideal for a wide range of healthcare environments and treatments procedures.

**TMBB** - Two Section Electric Height Elevation

**TMBM** - Single Section Electric Height Elevation

## STANDARD SPECIFICATION

- Safe working load and lifting capacity of 350kg (55 stone)
- Precision bearings on all pivot points
- Hand switch fitted as standard
- Power assisted head section with anti-trap design
- Excellent height range from 45cm (18") to 106cm (42")
- Hoist clearance along entire length of couch
- Gas spring assisted back rest to +85°
- Breathing hole and removable plug
- 50mm high density flame-retardant foam
- Antimicrobial contoured upholstery
- Fully compliant and CE marked

Safe working  
load and lifting  
capacity of 350kg  
(55 stone)

## OPTIONAL ACCESSORIES

- Motorised backrest
- Single foot switch operation for electric models
- Double foot switch operation for electric models
- Paper towel holder
- Single piece top section

## DIMENSIONS


- Length: 190cm
- Width: 101cm
- Safe Load: 350kg
- Height range: 45cm - 106cm

## WARRANTY

- 3 year warranty on moving components
- Lifetime warranty on frame

## STANDARDS

- EEC classification - MDD93/42/EEC Annex IX Class I
- Electrical standard - IEC601 Part 1 + 60601. P (EN)
- IP rating - IPx4
- Vinyl flame retardant to CRIB5 (BS7176 inc. foam)
- Foam to BS 5852: Part 2: 1982, Ignition Source 5 (CRIB5)



Width 101cm



124cm

65cm



Telephone: +44 (0) 1484 865586, Email: [info@titanmed.co.uk](mailto:info@titanmed.co.uk), Web: [www.titanmed.co.uk](http://www.titanmed.co.uk)  
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW  
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288

# BARIATRIC PODIATRY CHAIR



A very strong, robust and aesthetically pleasing podiatry couch designed to meet the increasing need to cater for the heavier patient with electric height and backrest functions as standard. Our bariatric podiatry couch features 55mm retractable positioning casters with an adjustable self-levelling foot making them ideal for a wide range of healthcare environments and treatment procedures.

**TMBP** - Bariatric podiatry electric height elevation

## STANDARD SPECIFICATION

- Excellent height range from 48cm (18") to 97cm (38")
- Split and extending leg section
- Safe working load and lifting capacity of 350kg (55 stone)
- Motorised back rest from -25° to +85°
- Stay level swing back arm rests for easy patient access (arms also removable)
- Purpose made precision bearings on all pivot points
- Electric height elevation
- The power assisted head section incorporates anti-trap design
- Hand switch fitted as standard
- Hoist clearance along entire length of couch
- 50mm high density flame-retardant foam
- Antimicrobial contoured upholstery
- Fully compliant and CE marked



## OPTIONAL ACCESSORIES


- Double foot switch
- Paper towel holder
- Breathing hole and plug

## WARRANTY

- 3 year warranty on moving components
- Lifetime warranty on frame

## DIMENSIONS

- Length: 190cm
- Width: 76cm
- Height range: 48cm - 97cm
- Safe Load: 350kg



Safe working  
load and lifting  
capacity of 350kg  
(55 stone)



## STANDARDS

- IP rating - IPx4
- EEC classification - MDD93/42/EEC Annex IX Class I
- Electrical standard - IEC601 Part 1 + 60601. P (EN)
- Vinyl flame retardant to CRIB5 (BS7176 inc. foam)
- Foam to BS 5852: Part 2: 1982, Ignition Source 5 (CRIB5)



Telephone: +44 (0) 1484 865586, Email: [info@titanmed.co.uk](mailto:info@titanmed.co.uk), Web: [www.titanmed.co.uk](http://www.titanmed.co.uk)  
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW  
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288

# NON-TILTING PODIATRY CHAIR



This entry level podiatry couch is a very strong, robust and aesthetically pleasing model making it ideal for a variety of healthcare environments and treatment procedures.

## TM3P - Double Footswitch Electric Podiatry Chair

### STANDARD SPECIFICATION

- Electric height adjustment and motorised backrest
- Leg section adjustable from  $-30^{\circ}$  to  $+30^{\circ}$ , seat tilt from  $0^{\circ}$  to  $+25^{\circ}$  & backrest adjustment from  $-15^{\circ}$  to  $+90^{\circ}$
- Stay level swing back arm rests for easy patient access (arms also removable)
- Leg length adjustment  $+15\text{cm}$  (6") on divided leg model
- Robust retractable castors for excellent mobility and safety
- Purpose made precision bearings on all pivot joints
- Footrest adjustment by gas strut mechanism
- Adjustable foot for uneven floors
- Couch goes horizontal for patient recovery in emergency situations

### OPTIONAL ACCESSORIES

- Single Leg Section
- Breathing hole and plug
- Paper roll holder
- Double thickness padded upholstery
- Matching operators chair or stool

### WARRANTY


- 3 year warranty on moving components
- Lifetime warranty on frame
- Fully compliant and CE marked

### STANDARDS


- EEC classification - MDD93/42/EEC Annex IX Class I
- Electrical standard - IEC601 Part 1 + 60601. P (EN)
- Vinyl flame retardant to CRIB5 (BS7176 inc. foam)
- Foam to BS 5852: Part 2: 1982, Ignition Source 5 (CRIB5)
- IP rating - IPx4

### DIMENSIONS

- Length: 174cm x Width: 63cm
- Height range: 48cm - 97cm
- Safe Load: 250kg



Height Range 48cm - 97cm



Width 63cm



Safe working  
load and lifting  
capacity of 250kg  
(39 stone)



Telephone: +44 (0) 1484 865586, Email: [info@titanmed.co.uk](mailto:info@titanmed.co.uk), Web: [www.titanmed.co.uk](http://www.titanmed.co.uk)  
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW  
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288

# TILTING PODIATRY CHAIR



Being aesthetically pleasing, very strong and robust, this triple motor chair is our most popular podiatry chair and has all features of the TM3P model but with the added benefit of Trendelenburg seat tilt and triple footswitch included as standard. This greatly enhances patient comfort and allows the foot end to rise up to a full 65" (165cm) for unrivalled access whilst doing heelwork or for operating when standing up.

## TM3PT - Triple Footswitch Electric Tilting Podiatry Chair

### STANDARD SPECIFICATION

- Motorised backrest and triple footswitch included as standard
- Electric height, motorised backrest and trendelenburg tilt adjustment
- Leg section adjustable from -30° to +30°, seat tilt from 0° to +25° & backrest adjustment from -15° to +90°
- Stay level swing back arm rests for easy patient access (arms also removable)
- Leg length adjustment +15cm (6") on divided leg model
- Purpose made precision bearings on all pivot joints
- Couch goes horizontal for patient recovery in emergency situations
- Robust retractable castors for excellent mobility and safety
- Footrest adjustment by gas strut mechanism
- Adjustable foot for uneven floors

Safe working  
load and lifting  
capacity of 250kg  
(39 stone)

### OPTIONAL ACCESSORIES

- Single Leg Section
- Breathing hole and plug
- Paper roll holder
- Double thickness padded upholstery
- Matching operators chair or stool

### STANDARDS



- EEC classification - MDD93/42/EEC Annex IX Class I
- Electrical standard - IEC601 Part 1 + 60601. P (EN)
- Vinyl flame retardant to CRIB5 (BS7176 inc. foam)
- Foam to BS 5852: Part 2: 1982, Ignition Source 5 (CRIB5)
- IP rating - IPx4

### WARRANTY

- 3 year warranty on moving components
- Lifetime warranty on frame
- Fully compliant & CE marked

### DIMENSIONS

- Length: 174cm
- Width: 63cm
- Safe Load: 250kg
- Height range: 48cm - 97cm



Motorised  
backrest & triple  
footswitch as  
standard



Telephone: +44 (0) 1484 865586, Email: [info@titanmed.co.uk](mailto:info@titanmed.co.uk), Web: [www.titanmed.co.uk](http://www.titanmed.co.uk)  
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW  
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288

# DROP END PODIATRY CHAIR



This entry level single motor podiatry treatment chair comes with full 90° leg drop and a shorter seat section allowing unrivalled patient comfort and access. This model features smooth gas strut operated back rest and divided leg rests.

## TMDEP - Drop End Electric Podiatry Chair

### STANDARD SPECIFICATION

- Electric height adjustment
- Leg section adjustable from -90° to +5° & backrest adjustment from -15° to +90°
- Stay level swing back arm rests for easy patient access (arms also removable)
- Robust retractable castors for excellent mobility and safety
- Purpose made precision bearings on all pivot joints
- Footrest adjustment by gas strut mechanism
- Adjustable foot for uneven floors
- Couch goes horizontal for patient recovery in emergency situations

### OPTIONAL ACCESSORIES

- Motorised backrest
- Single leg section
- Breathing hole and plug
- Paper roll holder
- Double thickness padded upholstery
- Matching operators chair or stool

Safe working  
load and lifting  
capacity of 250kg  
(39 stone)

### STANDARDS


- EEC classification - MDD93/42/EEC Annex IX Class I
- IP rating - IPx4
- Electrical standard - IEC601 Part I + 60601. P (EN)
- Vinyl flame retardant to CRIB5 (BS7176 inc. foam)
- Foam to BS 5852: Part 2: 1982, Ignition Source 5 (CRIB5)

### WARRANTY

- 3 year warranty on moving components
- Lifetime warranty on frame
- Fully compliant and CE marked

### DIMENSIONS

- Length: 174cm
- Width: 62cm
- Height range: 48cm - 97cm
- Safe Load: 250kg



Telephone: +44 (0) 1484 865586, Email: [info@titanmed.co.uk](mailto:info@titanmed.co.uk), Web: [www.titanmed.co.uk](http://www.titanmed.co.uk)  
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW  
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288

# 2 SECTION BARIATRIC



A very strong, robust and aesthetically pleasing treatment couch designed to meet the increasing need to cater for the heavier patient with electric height and backrest functions as standard. Our two section bariatric couch features 55mm retractable positioning casters with an adjustable self-levelling foot making them ideal for a wide range of healthcare environments and treatment procedures.

**TMB2** - Electric height elevation

## STANDARD SPECIFICATION

- Safe working load and lifting capacity of 350kg (55 stone)
- Excellent height range from 43cm (17") to 97cm (38")
- Electric height elevation
- Antimicrobial contoured upholstery
- Fully compliant and CE marked
- 50mm high density flame-retardant foam
- The power assisted head section incorporates anti-trap design
- Hand switch fitted as standard
- Hoist clearance along entire length of couch
- Purpose made precision bearings on all pivot points
- Motorised back rest from -25° to +85°

## OPTIONAL ACCESSORIES

- Double foot switch operation
- Paper towel holder
- Phlebotomy arm
- Breathing hole and plug
- Cot sides (powder coated, chrome plated or padded)

## DIMENSIONS


- Length: 190cm
- Width: 76cm
- Safe Load: 350kg
- Height range: 43cm - 97cm

## WARRANTY

- 3 year warranty on moving components
- Lifetime warranty on frame

## STANDARDS

- EEC classification - MDD93/42/EEC Annex IX Class I
- Electrical standard - IEC601 Part 1 + 60601. P (EN)
- IP rating - IPx4
- Vinyl flame retardant to CRIB5 (BS7176 inc. foam)
- Foam to BS 5852: Part 2: 1982, Ignition Source 5 (CRIB5)



Telephone: +44 (0) 1484 865586, Email: [info@titanmed.co.uk](mailto:info@titanmed.co.uk), Web: [www.titanmed.co.uk](http://www.titanmed.co.uk)  
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW  
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288

# 3 SECTION BARIATRIC



A very strong, robust and aesthetically pleasing treatment couch designed to meet the increasing need to cater for the heavier patient with electric height and backrest functions as standard. Our three section bariatric couch features 55mm retractable positioning casters with an adjustable self-levelling foot making them ideal for a wide range of healthcare environments and treatment procedures.

**TMB3** - Electric height elevation

## STANDARD SPECIFICATION

- Safe working load and lifting capacity of 350kg (55 stone)
- Electric height elevation
- Motorised back rest from -25° to +85°
- Purpose made precision bearings on all pivot points
- Hoist clearance along entire length of couch
- Excellent height range from 48cm (18") to 97cm (38")
- The power assisted head section incorporates anti-trap design
- Hand switch fitted as standard
- 50mm high density flame-retardant foam
- Antimicrobial contoured upholstery
- Fully compliant and CE marked



## OPTIONAL ACCESSORIES

- Double foot switch
- Paper towel holder
- Phlebotomy arm
- Breathing hole and plug
- Cot sides (powder coated, chrome plated or padded)

Safe working  
load and lifting  
capacity of 350kg  
(55 stone)

## DIMENSIONS

- Width: 76cm
- Height range: 48cm - 97cm
- Length: 190cm
- Safe Load: 350kg



## WARRANTY

- 3 year warranty on moving components
- Lifetime warranty on frame

## STANDARDS

- EEC classification - MDD93/42/EEC Annex IX Class I
- IP rating - IPx4
- Electrical standard - IEC601 Part I + 60601. P (EN)
- Vinyl flame retardant to CRIB5 (BS7176 inc. foam)
- Foam to BS 5852: Part 2: 1982, Ignition Source 5 (CRIB5)



Telephone: +44 (0) 1484 865586, Email: [info@titanmed.co.uk](mailto:info@titanmed.co.uk), Web: [www.titanmed.co.uk](http://www.titanmed.co.uk)  
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW  
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288

# GYNAECOLOGY COUCH



A treatment couch designed for general gynaecology procedures featuring electric adjustment of height, a fixed 5° tilted seat and manual backrest. All models feature a crescent shaped seat section, stainless steel debris tray, arm loops and fully adjustable leg supports.

**TMGIE** - Gynaecology Couch - electric height elevation

## STANDARD SPECIFICATION

- Safe working load and lifting capacity of 250kg (39 stone)
- Electric height adjustment from 52cm (20") to 106cm (42")
- Extension piece with shaped pad to fit in crescent cut out
- Adjustable manual backrest +85°
- Compact base with retractable castors
- Antimicrobial contoured upholstery
- Precision bearings on all pivot points
- 50mm high density flame-retardant foam
- Hoist clearance along length of couch
- Adjustable foot for stability on uneven ground

Safe working  
load and lifting  
capacity of 250kg  
(39 stone)

## OPTIONAL ACCESSORIES

- Seat extension
- Paper roll holder

## DIMENSIONS


- Length: 125cm (with extension piece 163cm)
- Width: 65cm
- Safe Load: 250kg
- Height range: 52cm - 106cm

## WARRANTY

- 3 year warranty on moving components
- Lifetime warranty on frame
- Fully compliant and CE marked

## STANDARDS

- EEC classification - MDD93/42/EEC Annex IX Class I
- Electrical standard - IEC601 Part 1 + 60601. P (EN)
- IP rating - IPx4
- Vinyl flame retardant to CRIB5 (BS7176 inc. foam)
- Foam to BS 5852: Part 2: 1982, Ignition Source 5 (CRIB5)



Width 65cm



52cm

72cm

Height Range 52cm - 106cm



Telephone: +44 (0) 1484 865586, Email: [info@titanmed.co.uk](mailto:info@titanmed.co.uk), Web: [www.titanmed.co.uk](http://www.titanmed.co.uk)  
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW  
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288

# GYNAECOLOGY COUCH



A treatment couch designed for general gynaecology procedures featuring full electric adjustment of height and tilt. All models feature a crescent shaped seat section, stainless steel debris tray, arm loops and fully adjustable leg supports.

**TMG3E** - Gynaecology Couch - electric height elevation, tilt and backrest

## STANDARD SPECIFICATION

- Safe working load and lifting capacity of 250kg (39 stone)
- Electric height adjustment from 52cm (20") to 106cm (42")
- Adjustable electric backrest +85°
- Extension piece with shaped pad to fit in crescent cut out
- Compact base with retractable castors
- Adjustable foot for stability on uneven ground
- Antimicrobial contoured upholstery
- Precision bearings on all pivot points
- 50mm high density flame-retardant foam
- Hoist clearance along length of couch
- Triple foot switch operation

## OPTIONAL ACCESSORIES

- Seat extension
- Paper roll holder

## DIMENSIONS



- Length: 125cm (with extension piece 163cm)
- Width: 65cm
- Safe Load: 250kg
- Height range: 52cm - 106cm

## WARRANTY

- 3 year warranty on moving components
- Lifetime warranty on frame

## STANDARDS

- EEC classification - MDD93/42/EEC Annex IX Class I
- Electrical standard - IEC601 Part I + 60601. P (EN)
- IP rating - IPx4
- Vinyl flame retardant to CRIB5 (BS7176 inc. foam)
- Foam to BS 5852: Part 2: 1982, Ignition Source 5 (CRIB5)



Telephone: +44 (0) 1484 865586, Email: [info@titanmed.co.uk](mailto:info@titanmed.co.uk), Web: [www.titanmed.co.uk](http://www.titanmed.co.uk)  
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW  
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288

# PHLEBOTOMY ELECTRIC COUCH



With the foot section dropping to 90°, this versatile and compact drop end three section couch allows for optimum use in smaller treatment areas. A specialist couch designed to provide a comfortable and easy access seat position for blood sampling procedures. All models are fitted with a pair of multi-adjustable arm rests which can be easily removed allowing the chair to be used for a range of minor procedures.

**TMPB** - Phlebotomy Couch electric height elevation

## STANDARD SPECIFICATION

- Electric height adjustment
- Leg section adjustable from -90° to +5°
- Backrest adjustment from -15° to +90°
- Stay level swing back arm rests for easy patient access (arms also removable)
- Couch goes horizontal for patient recovery in emergency situations
- Robust retractable castors for excellent mobility and safety
- Purpose made precision bearings on all pivot joints
- Footrest adjustment by gas strut mechanism
- Adjustable foot for uneven floors

## OPTIONAL ACCESSORIES

- Motorised backrest
- Single leg section
- Breathing hole and plug
- Paper roll holder
- Double thickness padded upholstery

## STANDARDS

- EEC classification - MDD93/42/EEC Annex IX Class I
- IP rating - IPx4
- Electrical standard - IEC601 Part 1 + 60601. P (EN)
- Vinyl flame retardant to CRIB5 (BS7176 inc. foam)
- Foam to BS 5852: Part 2: 1982, Ignition Source 5 (CRIB5)

## WARRANTY


- 3 year warranty on moving components
- Lifetime warranty on frame
- Fully compliant and CE marked

## DIMENSIONS

- Length: 174cm
- Width: 62cm
- Height range: 48cm - 97cm
- Safe Load: 250kg



Safe working  
load and lifting  
capacity of 250kg  
(39 stone)



Telephone: +44 (0) 1484 865586, Email: [info@titanmed.co.uk](mailto:info@titanmed.co.uk), Web: [www.titanmed.co.uk](http://www.titanmed.co.uk)  
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW  
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288

# ULTRASOUND MEDICAL COUCH



An electric height and motorised three section ultrasound couch featuring a Trendelenburg tilt to allow for patient scanning. The couch features 55mm retractable castors with locking function as well as a 90° drop for patient access.

## TMUS - Ultrasound Medical Couch

### STANDARD SPECIFICATION

- Electric height, motorised backrest and trendelenburg tilt adjustment
- Leg section adjustable to -90°, seat tilt from 0° to +25° & backrest adjustment from -15° to +90°
- Stay level swing back arm rests for easy patient access (arms also removable)
- Leg length adjustment +15cm (6") on divided leg model
- Couch goes horizontal for patient recovery in emergency situations
- Purpose made precision bearings on all pivot joints
- Footrest adjustment by gas strut mechanism
- Adjustable foot for uneven floors
- Retractable castor wheels with locking function

### OPTIONAL ACCESSORIES

- Breathing hole and plug
- Paper roll holder
- Double thickness padded upholstery
- Matching operators chair or stool

### STANDARDS


- IP rating - IPx4
- EEC classification - MDD93/42/EEC Annex IX Class I
- Electrical standard - IEC601 Part I + 60601. P (EN)
- Vinyl flame retardant to CRIB5 (BS7176 inc. foam)
- Foam to BS 5852: Part 2: 1982, Ignition Source 5 (CRIB5)

### WARRANTY


- 3 year warranty on moving components
- Lifetime warranty on frame
- Fully compliant and CE marked

### DIMENSIONS

- Length: 174cm, Width: 63cm
- Height range: 45cm - 106cm
- Safe Load: 250kg



Height Range 45cm - 106cm



Width 63cm



Safe working  
load and lifting  
capacity of 250kg  
(39 stone)



Telephone: +44 (0) 1484 865586, Email: [info@titanmed.co.uk](mailto:info@titanmed.co.uk), Web: [www.titanmed.co.uk](http://www.titanmed.co.uk)  
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW  
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288

# DROP END MULTI - PURPOSE COUCH

**TITANMED**

This robust, versatile and compact drop end three section couch allows for optimum use of smaller treatment areas. With the foot section dropping to 90°, the couch is ideal for use in a wide range of healthcare and clinical environments.

**TMDEE** - Drop End Electric Multi Purpose Couch

**TMDEH** - Drop End Hydraulic Multi Purpose Couch

## STANDARD SPECIFICATION

- Choice of 2 models - hydraulic or electric height elevation
- Safe working load and lifting capacity of 250kg (39 stone)
- Gas spring assisted back rest from -30° to +85°
- Excellent height range from 48cm to 97cm
- Purpose made precision bearings on all pivot points
- Hoist clearance along entire length of couch
- Hand switch fitted as standard on electric models
- The power assisted head section incorporates anti-trap design
- 50mm high density flame-retardant foam
- Antimicrobial contoured upholstery
- Fully compliant and CE marked

## OPTIONAL ACCESSORIES

- Motorised backrest
- Single foot switch operation for electric models
- Double foot switch operation for electric models
- Paper towel holder
- Cot sides (powder coated, chrome plated or padded)
- Breathing hole and plug

## DIMENSIONS



- Length: 192cm (head 69cm, seat 50cm, foot 68cm)
- Width: 65cm
- Safe Load: 250kg
- Height range: 48cm - 97cm

## WARRANTY

- 3 year warranty on moving components
- Lifetime warranty on frame

## STANDARDS

- EEC classification - MDD93/42/EEC Annex IX Class I
- IP rating - IPx4
- Electrical standard - IEC601 Part I + 60601. P (EN)
- Vinyl flame retardant to CRIB5 (BS7176 inc. foam)
- Foam to BS 5852: Part 2: 1982, Ignition Source 5 (CRIB5)



Telephone: +44 (0) 1484 865586, Email: [info@titanmed.co.uk](mailto:info@titanmed.co.uk), Web: [www.titanmed.co.uk](http://www.titanmed.co.uk)  
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW  
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288

# FIXED HEIGHT COUCH



Our fixed height surgery couch is available with optional drawer units and mobile storage cabinet in four contemporary finishes. The couch combines a number of treatment and storage solutions with excellent comfort for the client.

Available in a range of 20 upholstery colours, these examination couches are ideal for use in a wide range of healthcare and clinical environments, including GP surgeries, hospitals, first aid rooms and health centres.

## TMFH - Fixed Height Couch

### STANDARD SPECIFICATION

- Multi-position self locking backrest
- 2 Layer/density foam for support & comfort
- Luxury 80mm deep contoured upholstery
- Anti-bacterial/anti-microbial vinyl
- Paper roll holder supplied as standard
- Durable silver epoxy coated frame
- 40mm diameter legs for stability
- 1 Adjustable foot for uneven floors
- Choice of integrated storage solutions
- Easy clean design and materials

### DIMENSIONS

- Length: 185cm
- Width: 60cm
- Height: 82cm
- Safe Load: 185kg

### WARRANTY

- 2 year warranty on moving components
- Lifetime warranty on frame
- Fully compliant and CE marked

### STANDARDS

- EEC classification - MDD93/42/EEC Annex IX Class I
- Vinyl flame retardant to CRIB5 (BS7176 inc. foam)
- Foam to BS 5852: Part 2: 1982, Ignition Source 5 (CRIB5)

### COLOURS

- Choice of 3 wood colour options
- Choice of 20 fabric colour options



Telephone: +44 (0) 1484 865586, Email: [info@titanmed.co.uk](mailto:info@titanmed.co.uk), Web: [www.titanmed.co.uk](http://www.titanmed.co.uk)  
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW  
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288

# APOLLO PATIENT TROLLEY



This Patient Trolley is suited to patient examination, transportation, treatment and recovery in a range of healthcare departments.

It has easily accessible pedals to adjust height and

Trendelenburg and a simple, user-friendly release lever; the backrest is also entirely adjustable to 90°.

A safe working load of 250kg comes as standard.

## TMPS3 - Apollo Patient Trolley

### STANDARD SPECIFICATION

- Nano ceramic, electrostatic epoxy polyester powder metal parts
- 2-piece, easy clean HPL lying surface
- Single button collapsible side rails
- Bumpers and holes for accessories in all corners
- Height adjustable IV pole
- Stainless steel pushing bar located at both ends of stretcher
- 360° turning swivel casters equipped with central and directional lock (150mm)
- Gas spring controlled Trendelenburg adjustment
- Gas spring controlled backrest adjustment
- Oxygen bottle holder
- 8cm thick mattress suitable to the lying surface

### OPTIONAL ACCESSORIES

- Defibrillator table
- X-ray cassette carrier

### WARRANTY

- 2 year warranty on moving components

### STANDARDS


- EEC classification - MDD93/42/EEC Annex IX Class I
- Vinyl flame retardant to CRIB5 (BS7176 inc. foam)

### DIMENSIONS

- Width (Overall): 84.5cm x Length (Overall): 208.5cm
- Height: 60cm
- Lying Surface: 193.5 x 60cm
- Backrest Angle: 90°
- Trendelenburg: 12°



Safe working  
load and lifting  
capacity of 250kg  
(39 stone)



Telephone: +44 (0) 1484 865586, Email: [info@titanmed.co.uk](mailto:info@titanmed.co.uk), Web: [www.titanmed.co.uk](http://www.titanmed.co.uk)  
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW  
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288

# APOLLO EMERGENCY PATIENT TROLLEY

**TITANMED**

The Apollo Emergency Patient Trolley is designed for patient treatments and transportation within a wide range of hospital environments. Featuring a 5th wheel steering system the trolley is highly manoeuvrable and has easily accessible pedals to adjust height and Trendelenburg system. With a simple, user-friendly release lever, the backrest is also entirely adjustable to 90°. X-ray cassette holder under the entire lying surface.

Safe working  
load and lifting  
capacity of 250kg  
(39 stone)

## TMP55 - Emergency Patient Trolley

### STANDARD SPECIFICATION

- 2-piece, easy clean HPL lying surface
- Single button collapsible side rails
- Stainless steel pushing bar at both ends of stretcher
- Aluminium telescopic columns
- ABS plastic bottom part covers
- Bumpers and holes for accessories in corners
- Nano ceramic, electrostatic epoxy polyester powder metal parts
- Backrest angle adjustment by gas spring
- Height adjustment, Trendelenburg & Reverse Trendelenburg positions achieved by gas springs operated by pedals
- 360° turning swivel casters equipped with central and directional lock (200mm)
- 5th caster 125mm (for easy movement of stretcher)
- Height adjustable IV pole
- Oxygen bottle holder
- Patient belt
- 8cm thick mattress compatible to the lying surface
- X-ray cassette holder under entire lying surface

### OPTIONAL ACCESSORIES



- Defibrillator table

### WARRANTY

- 2 year warranty on moving components

### DIMENSIONS

- Width (Overall): 84.5cm x Length (Overall): 208.5cm
- Height: Minimum 61cm, Maximum 89cm
- Lying Surface: 193.5 x 60cm
- Backrest Angle: 90°
- Trendelenburg: 10° & Reverse Trendelenburg: 10°



Telephone: +44 (0) 1484 865586, Email: [info@titanmed.co.uk](mailto:info@titanmed.co.uk), Web: [www.titanmed.co.uk](http://www.titanmed.co.uk)  
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW  
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288

# ELECTRIC TILT TABLE



Designed as a strong, robust platform suitable for a wide range of applications from rehabilitation through to examinations, this tilt table has a four button hand control permitting simultaneous control and adjustment of both height and tilt functions. Lowering to just wheelchair height and featuring electric height and tilt, large toe braked castors, adjustable harnesses and patient grab handles.

**TMTT** - Electric Tilt Table

## STANDARD SPECIFICATION

- Electric height & tilt operation
- Hardwearing anti-microbial vinyl
- Safe working load 160kg
- 100mm toe braked castors
- Hoist clearance 130mm
- Swing back patient grab handles
- Set of patient harnesses
- Standard patient foot board
- Fully washable and impermeable, to aid infection control
- Suitable for cleaning and disinfecting with dilute solutions of Hypochlorite (10,000ppm)

Safe working  
load and lifting  
capacity of 160kg  
(25 stone)

## STANDARDS

- IP Rating - IPx4
- EEC Classification - MDD93/42/EEC Annex IX Class I
- Electrical Standard - IEC601 Part 1 + 60601. P (EN)
- Vinyl Flame Retardant To CRIB 5 (BS7176 Including Foam)
- Foam To BS 5852: Part 2: 1982, Ignition Source 5 (CRIB 5)

## DIMENSIONS

- Height: 60cm - 105cm approx. (24" - 41") approx.
- Length: 185cm (73") approx.
- Width: 65cm (26") approx.
- Operational parameters: 85° of tilt capabilities from horizontal

## OPTIONAL ACCESSORIES

- Divided Foot Board
- Battery Back-up
- Angle Indicator
- D Cushion

## WARRANTY

- 3 year warranty on moving components
- Lifetime warranty on frame
- Fully compliant and CE marked



Telephone: +44 (0) 1484 865586, Email: [info@titanmed.co.uk](mailto:info@titanmed.co.uk), Web: [www.titanmed.co.uk](http://www.titanmed.co.uk)  
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW  
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288

# MEDICAL CHAIRS & STOOLS



TitanMed offer a selection of chairs and stools, with a wide range of upholstery colour options to match your couch and easy clean bases as standard. All chairs and stools are stable and easily manoeuvrable, with five independently moving castors. Height adjustment is made simple and continuous with a hydraulic compression spring. Some models offer adjustable backrests, whilst others come with fixed backrests, allowing you to always find your perfect sitting position.

## Titan Operators Chair

- Adjustable height, tilt and backrest
- 63cm 5 castor base for maximum stability
- PU Upholstery

### DIMENSIONS

- Seat Depth: 44cm (17") approx.
- Seat Width: 45cm (18") Approx.
- Height: 44 -58cm (17" - 23") Approx.
- Backrest Height: 41cm (61") Approx.
- Backrest Width: 44cm (17") Approx.
- Diameter: 63cm (25") Approx.

**Code: OPSC**

Full choice  
of colours to  
match your  
couch



## Titan Saddle Stool

- Adjustable height and tilt
- 5 smooth rolling wheels
- PU Upholstery

### DIMENSIONS

- Seat width: 44cm (17") approx.
- Height: 57-75cm min-max (22-30") approx.

**Code: SADDLE**

Strengthens  
the spine and  
improves  
posture



## Titan Backrest Saddle Stool

- Adjustable height, tilt and backrest
- 5 smooth rolling wheels
- PU Upholstery

### DIMENSIONS

- Seat width: 44cm (17") approx.
- Height: 57-75cm min-max (22-30") approx.
- Width Backrest: 44cm

**Code: SADDLEBR**

Can relieve  
aches and pain in  
arms, shoulders  
& back



Telephone: +44 (0) 1484 865586, Email: [info@titanmed.co.uk](mailto:info@titanmed.co.uk), Web: [www.titanmed.co.uk](http://www.titanmed.co.uk)  
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW  
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288