

MEDICAL CHAIRS & STOOLS



TitanMed offer a selection of chairs and stools, with a wide range of upholstery colour options to match your couch and easy clean bases as standard. All chairs and stools are stable and easily manoeuvrable, with five independently moving castors. Height adjustment is made simple and continuous with a hydraulic compression spring. Some models offer adjustable backrests, whilst others come with fixed backrests, allowing you to always find your perfect sitting position.

Titan Operators Chair

- Adjustable height, tilt and backrest
- 63cm 5 castor base for maximum stability
- PU Upholstery

DIMENSIONS

- Seat Depth: 44cm (17") approx.
- Seat Width: 45cm (18") Approx.
- Height: 44 -58cm (17" - 23") Approx.
- Backrest Height: 41cm (61") Approx.
- Backrest Width: 44cm (17") Approx.
- Diameter: 63cm (25") Approx.

Code: OPSC

Full choice
of colours to
match your
couch



Titan Saddle Stool

- Adjustable height and tilt
- 5 smooth rolling wheels
- PU Upholstery

DIMENSIONS

- Seat width: 44cm (17") approx.
- Height: 57-75cm min-max (22-30") approx.

Code: SADDLE

Strengthens
the spine and
improves
posture



Titan Backrest Saddle Stool

- Adjustable height, tilt and backrest
- 5 smooth rolling wheels
- PU Upholstery

DIMENSIONS

- Seat width: 44cm (17") approx.
- Height: 57-75cm min-max (22-30") approx.
- Width Backrest: 44cm

Code: SADDLEBR

Can relieve
aches and pain in
arms, shoulders
& back



Made in the UK



Telephone: +44 (0) 1484 865586, Email: info@titanmed.co.uk, Web: www.titanmed.co.uk
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288