

NON-TILTING PODIATRY CHAIR



This entry level podiatry couch is a very strong, robust and aesthetically pleasing model making it ideal for a variety of healthcare environments and treatment procedures.

TM3P - Double Footswitch Electric Podiatry Chair

STANDARD SPECIFICATION

- Electric height adjustment and motorised backrest
- Leg section adjustable from -30° to $+30^{\circ}$, seat tilt from 0 to $+25^{\circ}$ & backrest adjustment from -15° to $+90^{\circ}$
- Stay level swing back arm rests for easy patient access (arms also removable)
- Leg length adjustment $+15\text{cm}$ (6") on divided leg model
- Robust retractable castors for excellent mobility and safety
- Purpose made precision bearings on all pivot joints
- Footrest adjustment by gas strut mechanism
- Adjustable foot for uneven floors
- Couch goes horizontal for patient recovery in emergency situations

OPTIONAL ACCESSORIES

- Single Leg Section
- Breathing hole and plug
- Paper roll holder
- Double thickness padded upholstery
- Matching operators chair or stool

WARRANTY


- 3 year warranty on moving components
- Lifetime warranty on frame
- Fully compliant and CE marked

STANDARDS


- EEC classification - MDD93/42/EEC Annex IX Class I
- Electrical standard - IEC601 Part 1 + 60601. P (EN)
- Vinyl flame retardant to CRIB5 (BS7176 inc. foam)
- Foam to BS 5852: Part 2: 1982, Ignition Source 5 (CRIB5)
- IP rating - IPx4

DIMENSIONS

- Length: 174cm x Width: 63cm
- Height range: 48cm - 97cm
- Safe Load: 250kg



Height Range 48cm - 97cm



Width 63cm



Safe working
load and lifting
capacity of 250kg
(39 stone)



Telephone: +44 (0) 1484 865586, Email: info@titanmed.co.uk, Web: www.titanmed.co.uk
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288