

TILTING PODIATRY CHAIR



Being aesthetically pleasing, very strong and robust, this triple motor chair is our most popular podiatry chair and has all features of the TM3P model but with the added benefit of Trendelenburg seat tilt and triple footswitch included as standard. This greatly enhances patient comfort and allows the foot end to rise up to a full 65" (165cm) for unrivalled access whilst doing heelwork or for operating when standing up.

TM3PT - Triple Footswitch Electric Tilting Podiatry Chair

STANDARD SPECIFICATION

- Motorised backrest and triple footswitch included as standard
- Electric height, motorised backrest and trendelenburg tilt adjustment
- Leg section adjustable from -30° to +30°, seat tilt from 0° to +25° & backrest adjustment from -15° to +90°
- Stay level swing back arm rests for easy patient access (arms also removable)
- Leg length adjustment +15cm (6") on divided leg model
- Purpose made precision bearings on all pivot joints
- Couch goes horizontal for patient recovery in emergency situations
- Robust retractable castors for excellent mobility and safety
- Footrest adjustment by gas strut mechanism
- Adjustable foot for uneven floors

Safe working
load and lifting
capacity of 250kg
(39 stone)

OPTIONAL ACCESSORIES

- Single Leg Section
- Breathing hole and plug
- Paper roll holder
- Double thickness padded upholstery
- Matching operators chair or stool

STANDARDS


- EEC classification - MDD93/42/EEC Annex IX Class I
- Electrical standard - IEC601 Part 1 + 60601. P (EN)
- Vinyl flame retardant to CRIB5 (BS7176 inc. foam)
- Foam to BS 5852: Part 2: 1982, Ignition Source 5 (CRIB5)
- IP rating - IPx4

WARRANTY


- 3 year warranty on moving components
- Lifetime warranty on frame
- Fully compliant & CE marked

DIMENSIONS

- Length: 174cm
- Width: 63cm
- Safe Load: 250kg
- Height range: 48cm - 97cm



Height Range 48cm - 97cm



15cm 33cm 52cm 74cm

Motorised
backrest & triple
footswitch as
standard



Telephone: +44 (0) 1484 865586, Email: info@titanmed.co.uk, Web: www.titanmed.co.uk
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288