

NEUROLOGY COUCH



A very strong, robust and aesthetically pleasing treatment couch with breathing hole and removable plug fitted as standard. Our two section neurology couch features 55mm retractable positioning castors with locking function and an adjustable self-levelling foot making them ideal for a wide range of healthcare environments and treatments procedures.

TMBB - Two Section Electric Height Elevation

TMBM - Single Section Electric Height Elevation

STANDARD SPECIFICATION

- Safe working load and lifting capacity of 350kg (55 stone)
- Precision bearings on all pivot points
- Hand switch fitted as standard
- Power assisted head section with anti-trap design
- Excellent height range from 45cm (18") to 106cm (42")
- Hoist clearance along entire length of couch
- Gas spring assisted back rest to +85 degrees
- Breathing hole and removable plug
- 50mm high density flame-retardant foam
- Antimicrobial contoured upholstery
- Fully compliant and CE marked

Safe working
load and lifting
capacity of 350kg
(55 stone)

OPTIONAL ACCESSORIES

- Motorised backrest
- Single foot switch operation for electric models
- Double foot switch operation for electric models
- Paper towel holder
- Single piece top section

DIMENSIONS



- Length: 190cm
- Width: 101cm
- Safe Load: 350kg
- Height range: 45cm - 106cm

WARRANTY

- 3 year warranty on moving components
- Lifetime warranty on frame

STANDARDS

- EEC classification - MDD93/42/EEC Annex IX Class I
- Electrical standard - IEC601 Part 1 + 60601. P (EN)
- IP rating - IPx4
- Vinyl flame retardant to CRIB5 (BS7176 inc. foam)
- Foam to BS 5852: Part 2: 1982, Ignition Source 5 (CRIB5)



Telephone: +44 (0) 1484 865586, Email: info@titanmed.co.uk, Web: www.titanmed.co.uk
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288