

DROP END MULTI - PURPOSE COUCH

TITANMED

This robust, versatile and compact drop end three section couch allows for optimum use of smaller treatment areas. With the foot section dropping to 90°, the couch is ideal for use in a wide range of healthcare and clinical environments.

TMDEE - Drop End Electric Multi Purpose Couch

TMDEH - Drop End Hydraulic Multi Purpose Couch

STANDARD SPECIFICATION

- Choice of 2 models - hydraulic or electric height elevation
- Safe working load and lifting capacity of 250kg (39 stone)
- Gas spring assisted back rest from -30° to +85°
- Excellent height range from 48cm to 97cm
- Purpose made precision bearings on all pivot points
- Hoist clearance along entire length of couch
- Hand switch fitted as standard on electric models
- The power assisted head section incorporates anti-trap design
- 50mm high density flame-retardant foam
- Antimicrobial contoured upholstery
- Fully compliant and CE marked

OPTIONAL ACCESSORIES

- Motorised backrest
- Single foot switch operation for electric models
- Double foot switch operation for electric models
- Paper towel holder
- Cot sides (powder coated, chrome plated or padded)
- Breathing hole and plug

DIMENSIONS



- Length: 192cm (head 69cm, seat 50cm, foot 68cm)
- Width: 65cm
- Safe Load: 250kg
- Height range: 48cm - 97cm

WARRANTY

- 3 year warranty on moving components
- Lifetime warranty on frame

STANDARDS

- EEC classification - MDD93/42/EEC Annex IX Class I
- IP rating - IPx4
- Electrical standard - IEC601 Part I + 60601. P (EN)
- Vinyl flame retardant to CRIB5 (BS7176 inc. foam)
- Foam to BS 5852: Part 2: 1982, Ignition Source 5 (CRIB5)



Telephone: +44 (0) 1484 865586, Email: info@titanmed.co.uk, Web: www.titanmed.co.uk
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288