

# ELITE RANGE

## 2 SECTION



Our top-of-the-range two section Elite treatment couch is both strong and robust whilst being aesthetically pleasing. Featuring twin lift mechanism for added stability, the Elite range is fitted with Power-Packer hydraulic pumps and Linak motors as standard. Ideal for a wide range of healthcare environments and treatments procedures with 55mm retractable positioning castors and an adjustable self-levelling foot.

**TME2H** - Hydraulic height elevation    **TME2E** - Electric height elevation

### STANDARD SPECIFICATION

- Safe working load and lifting capacity of 275kg (43 stone)
- Hydraulic or electric height elevation
- Excellent height range from 48cm to 102cm
- 40mm high density flame-retardant foam
- Antimicrobial contoured upholstery
- Fully compliant and CE marked
- Power assisted head section incorporates anti-trap design
- Hand switch fitted as standard on electric models
- Hoist clearance 135mm along entire length of couch
- Purpose made precision bearings on all pivot points
- Gas spring assisted back rest from -20° to +80°

### OPTIONAL ACCESSORIES

- Breathing hole with removable plug
- Single foot switch operation for electric models
- Double foam
- Paper towel holder
- Battery backup
- 75cm width

### DIMENSIONS

- Width: 63cm
- Length: 195cm (head 71cm, seat 124cm)
- Height range: 48cm - 102cm
- Safe Load: 275kg

### WARRANTY

- 3 year warranty on moving components
- Lifetime warranty on frame

### STANDARDS

- EEC classification - MDD93/42/EEC Annex IX Class I
- IP rating - IPx4
- Electrical standard - IEC601 Part I + 60601.P (EN)
- Vinyl flame retardant to CRIB5 (BS7176 inc. foam)
- Foam to BS5852: Part 2: 1982, Ignition Source 5 (CRIB5)
- Suitable for high traffic contract areas rated to 100,000 martindale rubs



Safe working  
load and lifting  
capacity of 275kg  
(43 stone)

Power-Packer  
hydraulic pumps  
& Linak motors  
as standard




Width 63cm




124cm



71cm



Height Range 48cm - 102cm



Telephone: +44 (0) 1484 865586, Email: [info@titanmed.co.uk](mailto:info@titanmed.co.uk), Web: [www.titanmed.co.uk](http://www.titanmed.co.uk)  
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW  
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288