

ELITE RANGE

3 SECTION



Our top-of-the-range three section Elite treatment couch is both strong and robust whilst being aesthetically pleasing. Featuring twin lift mechanism for added stability, the Elite range is fitted with Power-Packer hydraulic pumps and Linak motors as standard. Ideal for a wide range of healthcare environments and treatments procedures with 55mm retractable positioning castors and an adjustable self-levelling foot.

TME3H - Hydraulic height elevation **TME3E** - Electric height elevation

STANDARD SPECIFICATION

- Safe working load and lifting capacity of 275kg (43 stone)
- Gas spring assisted back rest from -25° to $+80^{\circ}$
- Hydraulic or electric height elevation
- Power assisted head section incorporates anti-trap design
- Excellent height range from 48cm to 102cm
- Purpose made precision bearings on all pivot points
- Hoist clearance 135mm along entire length of couch
- Hand switch fitted as standard on electric models
- 40mm high density flame-retardant foam
- Antimicrobial contoured upholstery
- Fully compliant and CE marked

OPTIONAL ACCESSORIES

- Breathing hole
- Single foot switch operation for electric models
- Double foam
- Paper towel holder
- Battery backup
- 75cm width

DIMENSIONS

- Width: 63cm
- Length: 195cm (head 71cm, seat 50cm, foot 71cm)
- Height range: 48cm - 102cm
- Safe Load: 275kg

WARRANTY

- 3 year warranty on moving components
- Lifetime warranty on frame

STANDARDS




- EEC classification - MDD93/42/EEC Annex IX Class I
- IP rating - IPx4
- Electrical standard - IEC601 Part I + 60601.P (EN)
- Vinyl flame retardant to CRIB5 (BS7176 inc. foam)
- Foam to BS5852: Part 2: 1982, Ignition Source 5 (CRIB5)
- Suitable for high traffic contract areas rated to 100,000 martindale rubs



Safe working load and lifting capacity of 275kg (43 stone)



Power-Packer hydraulic pumps & Linak motors as standard



Telephone: +44 (0) 1484 865586, Email: info@titanmed.co.uk, Web: www.titanmed.co.uk
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288