

FIXED HEIGHT COUCH



Our fixed height surgery couch is available with optional drawer units and mobile storage cabinet in four contemporary finishes. The couch combines a number of treatment and storage solutions with excellent comfort for the client.

Available in a range of 20 upholstery colours, these examination couches are ideal for use in a wide range of healthcare and clinical environments, including GP surgeries, hospitals, first aid rooms and health centres.

TMFH - Fixed Height Couch

STANDARD SPECIFICATION

- Multi-position self locking backrest
- 2 Layer/density foam for support & comfort
- Luxury 80mm deep contoured upholstery
- Anti-bacterial/anti-microbial vinyl
- Paper roll holder supplied as standard
- Durable silver epoxy coated frame
- 40mm diameter legs for stability
- 1 Adjustable foot for uneven floors
- Choice of integrated storage solutions
- Easy clean design and materials

DIMENSIONS

- Length: 185cm
- Width: 60cm
- Height: 82cm
- Safe Load: 185kg

WARRANTY

- 2 year warranty on moving components
- Lifetime warranty on frame
- Fully compliant and CE marked

STANDARDS

- EEC classification - MDD93/42/EEC Annex IX Class I
- Vinyl flame retardant to CRIB5 (BS7176 inc. foam)
- Foam to BS 5852: Part 2: 1982, Ignition Source 5 (CRIB5)

COLOURS

- Choice of 3 wood colour options
- Choice of 20 fabric colour options



Telephone: +44 (0) 1484 865586, Email: info@titanmed.co.uk, Web: www.titanmed.co.uk
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288