

# GYNAECOLOGY COUCH



A treatment couch designed for general gynaecology procedures featuring full electric adjustment of height and tilt. All models feature a crescent shaped seat section, stainless steel debris tray, arm loops and fully adjustable leg supports.

**TMG3E** - Gynaecology Couch - electric height elevation, tilt and backrest

## STANDARD SPECIFICATION

- Safe working load and lifting capacity of 250kg (39 stone)
- Electric height adjustment from 52cm (20") to 106cm (42")
- Adjustable electric backrest +85°
- Extension piece with shaped pad to fit in crescent cut out
- Compact base with retractable castors
- Adjustable foot for stability on uneven ground
- Antimicrobial contoured upholstery
- Precision bearings on all pivot points
- 50mm high density flame-retardant foam
- Hoist clearance along length of couch
- Triple foot switch operation

## OPTIONAL ACCESSORIES

- Seat extension
- Paper roll holder

## DIMENSIONS



- Length: 125cm (with extension piece 163cm)
- Width: 65cm
- Safe Load: 250kg
- Height range: 52cm - 106cm

## WARRANTY

- 3 year warranty on moving components
- Lifetime warranty on frame

## STANDARDS

- EEC classification - MDD93/42/EEC Annex IX Class I
- Electrical standard - IEC601 Part I + 60601. P (EN)
- IP rating - IPx4
- Vinyl flame retardant to CRIB5 (BS7176 inc. foam)
- Foam to BS 5852: Part 2: 1982, Ignition Source 5 (CRIB5)



Telephone: +44 (0) 1484 865586, Email: [info@titanmed.co.uk](mailto:info@titanmed.co.uk), Web: [www.titanmed.co.uk](http://www.titanmed.co.uk)  
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW  
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288