

PHLEBOTOMY ELECTRIC COUCH



With the foot section dropping to 90°, this versatile and compact drop end three section couch allows for optimum use in smaller treatment areas. A specialist couch designed to provide a comfortable and easy access seat position for blood sampling procedures. All models are fitted with a pair of multi-adjustable arm rests which can be easily removed allowing the chair to be used for a range of minor procedures.

TMPB - Phlebotomy Couch electric height elevation

STANDARD SPECIFICATION

- Electric height adjustment
- Leg section adjustable from -90° to +5°
- Backrest adjustment from -15° to +90°
- Stay level swing back arm rests for easy patient access (arms also removable)
- Couch goes horizontal for patient recovery in emergency situations
- Robust retractable castors for excellent mobility and safety
- Purpose made precision bearings on all pivot joints
- Footrest adjustment by gas strut mechanism
- Adjustable foot for uneven floors

OPTIONAL ACCESSORIES

- Motorised backrest
- Single leg section
- Breathing hole and plug
- Paper roll holder
- Double thickness padded upholstery

STANDARDS


- EEC classification - MDD93/42/EEC Annex IX Class I
- IP rating - IPx4
- Electrical standard - IEC601 Part I + 60601. P (EN)
- Vinyl flame retardant to CRIB5 (BS7176 inc. foam)
- Foam to BS 5852: Part 2: 1982, Ignition Source 5 (CRIB5)

WARRANTY

- 3 year warranty on moving components
- Lifetime warranty on frame
- Fully compliant and CE marked

DIMENSIONS

- Length: 174cm
- Width: 62cm
- Height range: 48cm - 97cm
- Safe Load: 250kg



Telephone: +44 (0) 1484 865586, Email: info@titanmed.co.uk, Web: www.titanmed.co.uk
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288