

APOLLO PATIENT TROLLEY



This Patient Trolley is suited to patient examination, transportation, treatment and recovery in a range of healthcare departments.

It has easily accessible pedals to adjust height and

Trendelenburg and a simple, user-friendly release lever; the backrest is also entirely adjustable to 90°.

A safe working load of 250kg comes as standard.

TMPS3 - Apollo Patient Trolley

STANDARD SPECIFICATION

- Nano ceramic, electrostatic epoxy polyester powder metal parts
- 2-piece, easy clean HPL lying surface
- Single button collapsible side rails
- Bumpers and holes for accessories in all corners
- Height adjustable IV pole
- Stainless steel pushing bar located at both ends of stretcher
- 360° turning swivel casters equipped with central and directional lock (150mm)
- Gas spring controlled Trendelenburg adjustment
- Gas spring controlled backrest adjustment
- Oxygen bottle holder
- 8cm thick mattress suitable to the lying surface

OPTIONAL ACCESSORIES

- Defibrillator table
- X-ray cassette carrier

WARRANTY

- 2 year warranty on moving components

STANDARDS



- EEC classification - MDD93/42/EEC Annex IX Class I
- Vinyl flame retardant to CRIB5 (BS7176 inc. foam)

DIMENSIONS

- Width (Overall): 84.5cm x Length (Overall): 208.5cm
- Height: 60cm
- Lying Surface: 193.5 x 60cm
- Backrest Angle: 90°
- Trendelenburg: 12°



Safe working
load and lifting
capacity of 250kg
(39 stone)



Telephone: +44 (0) 1484 865586, Email: info@titanmed.co.uk, Web: www.titanmed.co.uk
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288