

ELECTRIC TILT TABLE



Designed as a strong, robust platform suitable for a wide range of applications from rehabilitation through to examinations, this tilt table has a four button hand control permitting simultaneous control and adjustment of both height and tilt functions. Lowering to just wheelchair height and featuring electric height and tilt, large toe braked castors, adjustable harnesses and patient grab handles.

TMTT - Electric Tilt Table

STANDARD SPECIFICATION

- Electric height & tilt operation
- Hardwearing anti-microbial vinyl
- Safe working load 160kg
- 100mm toe braked castors
- Hoist clearance 130mm
- Swing back patient grab handles
- Set of patient harnesses
- Standard patient foot board
- Fully washable and impermeable, to aid infection control
- Suitable for cleaning and disinfecting with dilute solutions of Hypochlorite (10,000ppm)

Safe working
load and lifting
capacity of 160kg
(25 stone)

STANDARDS

- IP Rating - IPx4
- EEC Classification - MDD93/42/EEC Annex IX Class I
- Electrical Standard - IEC601 Part 1 + 60601. P (EN)
- Vinyl Flame Retardant To CRIB 5 (BS7176 Including Foam)
- Foam To BS 5852: Part 2: 1982, Ignition Source 5 (CRIB 5)

DIMENSIONS

- Height: 60cm - 105cm approx. (24" - 41") approx.
- Length: 185cm (73") approx.
- Width: 65cm (26") approx.
- Operational parameters: 85° of tilt capabilities from horizontal

OPTIONAL ACCESSORIES

- Divided Foot Board
- Battery Back-up
- Angle Indicator
- D Cushion

WARRANTY

- 3 year warranty on moving components
- Lifetime warranty on frame
- Fully compliant and CE marked



Telephone: +44 (0) 1484 865586, Email: info@titanmed.co.uk, Web: www.titanmed.co.uk
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288