

# 2 SECTION ELECTRIC & MANUAL



A very strong, robust and aesthetically pleasing treatment couch available with optional breathing hole and removable plug. Our two section couches feature 55mm retractable positioning castors with an adjustable self-levelling foot making them ideal for a wide range of healthcare environments and treatments procedures.

**TMU2H** - Hydraulic height elevation    **TMU2E** - Electric height elevation

## STANDARD SPECIFICATION

- Safe working load and lifting capacity of 250kg (39 stone)
- Choice of 2 models - hydraulic or electric height elevation
- Excellent height range from 48cm to 97cm
- 50mm high density flame-retardant foam
- Antimicrobial contoured upholstery
- Fully compliant and CE marked
- The power assisted head section incorporates anti-trap design
- Hand switch fitted as standard on electric models
- Hoist clearance along entire length of couch
- Purpose made precision bearings on all pivot points
- Gas spring assisted back rest from -30° to +85°

Safe working  
load and lifting  
capacity of 250kg  
(39 stone)



## OPTIONAL ACCESSORIES

- Motorised backrest
- Single foot switch operation for electric models
- Double foot switch operation for electric models
- Paper towel holder
- Cot sides (powder coated, chrome plated or padded)

## DIMENSIONS

- Length: 192cm (head 69cm, seat 122cm)
- Width: 65cm
- Safe Load: 250kg
- Height range: 48cm - 97cm

## WARRANTY

- 3 year warranty on moving components
- Lifetime warranty on frame

## STANDARDS


- EEC classification - MDD93/42/EEC Annex IX Class I
- Electrical standard - IEC601 Part 1 + 60601. P (EN)
- IP rating - IPx4
- Vinyl flame retardant to CRIB5 (BS7176 inc. foam)
- Foam to BS 5852: Part 2: 1982, Ignition Source 5 (CRIB5)



122cm



69cm



Height Range 48cm - 97cm



Telephone: +44 (0) 1484 865586, Email: [info@titanmed.co.uk](mailto:info@titanmed.co.uk), Web: [www.titanmed.co.uk](http://www.titanmed.co.uk)  
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW  
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288